

The art of seeing
what others can't.



Handcrafted Loupes, Lights & Camera

admetec
We've got your back



We've got your back.

Admetec handcrafted, custom-designed loupes, lights, and Camera

We specialize in creating solutions that are exceptionally lightweight, convenient to use, and ergonomic - enabling you to focus on being effective while we focus on your comfort.

Our handcrafted loupes, lights, and Camera are custom designed for YOU. To do so, we take into consideration your unique facial structure, working distance, style, and personal preferences. User-centered product design combined with high-quality materials and components is the foundation of Admetec excellence.

Seeing inside a patient's mouth is hard.
Sharing what you see is even harder.

Why Choose Flamingo

- ✓ **Effortless documentation** – Record Full HD video hands free, directly from your line of sight.
- ✓ **Powerful teaching tool** – Stream live to students or colleagues for real time learning.
- ✓ **Patient trust** – Show patients exactly what you see and strengthen informed consent.

Designed for Daily Practice

Dentists & Hygienists – Case documentation, social media content, patient education

Surgeons – Guided procedures, OR broadcasting, medicolegal recording

Educators & KOLs – Live demonstrations, course creation, conference streaming



Flamingo™ features

Feature	Benefit
Ultra lightweight (19 g)	Only 19g! - All day comfort-lighter than most loupe lights
Integrated 60 k lux LED	Powerful, shadow free illumination
Full HD 1080 p 30 fps camera	Crisp, stable footage for clinical documentation.
Remote control	Full control of the camera during operation, including the ability to take still images while recording and freeze the image for learning
Single hybrid cable	Truly wireless
Local real time Wi Fi broadcast	Share instantly on any PC or Mac-no internet required
AI Editing	AI Cleanup removes all unnecessary parts from the movie with one click
Swappable lenses (25 mm / 21 mm)	Select a narrow (25mm) or wide (21mm) field of view
Cloud storage	No need to rely on local memory-all videos are safely backed up to your account, under the user control
Up to 4 h camera + light or 8 h light only	Two swappable batteries keep you filming all day
Super Connector	360° rotation and vertical video mode (9:16 format)
Voice recording	Capture synchronized audio directly with Flamingo
Flamingo Software	Snapshots, essential editing, and free updates included
Compatible with most loupe brands	Adapters available for most of loupe brands
Data security	Fully compliant with HIPAA and GDPR requirements for patient privacy and data protection





We've got **your back.**

Admetec designs, develops, manufactures, and distributes handcrafted, custom-fitted magnifying loupes and LED headlamps as well as an integrated camera and headlight. We specialize in creating products that are exceptionally lightweight, convenient to use, and ergonomic – enabling you to focus on being effective while we focus on your comfort.

IM4000020 rev.B

